

Features

- · By using of button and a marking pen, 2 counter mode are possible.
- · This unit main power is 220V.
- · Easy to read 4 digits LED digital display.
- · Buzzer sounds on every counting.
- · A marking pen is usable for user's convenience.
- · Wider view of 125 mm wide lens.
- · Bacteria can be counted by picking with stainless needle point pen.

Specifications

Made in Korea

Spec. / Model	HYC-560
Display type	4 Digital LED Display
Petri dish holder	Accepts 100mm dish (removable) 150mm (fixed)
Zero rest button	Yes
Illumination	Fluorescent lamp
Dimensions (WxDxH,mm)	310x290x100
Mode of Counting	Manual push button counting
	Needle-point probing
	Marking pen probing
Electric Requirements	220VAC. 50/60Hz

Standard accessories	
Magnifier 5" 1ea	Stainless Needle point pen 1ea
Power cord 1ea	Making pen probe 1ea




Digital Colony Counter

www.hysc.co.kr

코로니카운터

Distinction



Structure

